

Catalogue 2025

primediagnostics

Droevendaalsesteeg 1
P.O. Box 16
6700 AA Wageningen
The Netherlands

WhatsApp: +31 317 480 613
E-mail : primediagnos@wur.nl
Web page: www.primediagnos.com
Web shop: shop.wur.nl/primediagnos



Table of content

Antisera Product information	1
General information on antisera.	1
Storage information and stability.	1
Information concerning the quality of the antisera.	1
Antisera to plant-pathogenic bacteria	1
Virus inocula	1
Orders and invoicing	2
Contact persons for information.	2
Abbreviations	2
Antisera Virus	3
Antisera Bacteria and Oomycetes	4
Prime Diagnostics Terms of Delivery	5
Workflow overview	6
Schematic overview of the ELISA procedure	6
Schematic overview of the Luminex xMAP procedure	7
Schematic overview of the TaqMan procedure	8
Schematic overview of the Lateral Flow Device procedure.	9
ELISA and Luminex xMAP comparison.	10
Appendix 1. Accessing webshop and creating an account	11

General information on antisera

- Antisera that are not on the list can potentially be produced on request. Please inquire.
- Custom-made conjugates with fluorescent labels or AP can be produced upon request.
- Protocols are available at www.primediagnosics.com.
- The AP-conjugate and coating should be diluted 1000x.
- The coating and alkaline phosphatase (AP) conjugates are suitable for use in DAS-ELISA only.

Storage information and stability

- Coatings and AP-conjugates should preferably be stored at 4 °C. At 4 °C antisera are stable for at least two years.
- Controls and IgG should be stored in aliquots at -20 °C. Avoid repeated freezing and thawing. When treated properly these reagents will remain stable for at least two years.
- For long time storage (periods longer than one year) store at -20 °C, prepare aliquots if necessary.
- When storing at -20 °C avoid repeated freezing and thawing, as this can influence the product quality considerably.
- Lyophilized positive control (VPC) contain lyophilized leaf material. The VPC remains stable for 1 year when stored at 4 °C. A prepared VPC can be stored in aliquots of 0.2 mL and is stable for 1 month at -20 °C and at 4 °C for 1 day.
- Lateral flow devices, packed in their original unopened (sealed) envelopes can be stored at room temperature. They will remain stable until expiration date. Lateral flow devices, packed in their original opened envelopes will remain stable for two weeks.
- Luminex reagents should be stored at 4°C, unless indicated different on the box. They will remain stable for one year.

Information concerning the quality of the antisera

- Antisera are biological products and differences in reactivity between batches are apparent. Prime Diagnostics exercises strict quality control and guarantees the reactivity as indicated on the batch quality sheet (if applicable).
- Negative controls are optimized for use in different crops and are not suitable for determination of the selectivity of the assay.
- Positive controls are qualitative and cannot be used for quantification!
- The criteria for reactivity of ELISA reagents for bacterial detection are an OD₄₀₅ reading of the 10 times diluted positive control (end concentration of 10⁷ cells/mL) of OD₄₀₅ > 1.0 after 30 minutes of substrate incubation.
- The criterion for reactivity of ELISA reagents for virus detection is an OD₄₀₅ reading of the 10 times diluted positive control of at least 0.5 OD after 30 minutes of substrate incubation. The value of the OD₄₀₅ may however vary depending on the ELISA reader used.
- Within our quality control procedure, our antisera are tested for lack of reactivity against a variety of pathogens and/or plant compounds known to occur in combination with the given pathogen. However, within the context of any given test, the possibility of false-positive results with unrelated organisms or plant or matrix substances should always be considered.
- Within our quality control procedure, our antisera are tested for specificity against a wide variety of currently known and recognized pathogenic isolates, strains and patho-types of the subsequent pathogens. However, new pathogenic or non-pathogenic strains and isolates of a given pathogen may emerge over time. Prime Diagnostics cannot guarantee that such a new strain or isolate will be detected with equal efficiency.

Antisera to plant-pathogenic bacteria

- Antisera labeled with different conjugates are available on request. Prices can be supplied upon request.
- For all bacteria non-infectious positive controls are available for use in DAS-ELISA or IIF. These controls have a concentration of 10⁸ cells/mL.
- Negative controls consisting of non-infectious bacteria or plant extract are available on request.
- The dilution of the IgG for IIF is dependent on the assay used and laboratory circumstances. It is recommended to determine the working dilution of the reagents under on-site conditions.
- The IgG is suitable for use in indirect immunofluorescence (IIF) detection.

Virus inocula

- Lyophilized virus particles, dried plant material (to be reconstituted in buffer) or purified virus for inoculating purposes can be obtained on request. The amounts are generally sufficient for small-scale screening experiments.

Orders and invoicing

WhatsApp : +31 317 480 613
 e-mail : primediagnosics@wur.nl
 web shop : shop.wur.nl/primediagnosics
 web page : www.primediagnosics.com

Contact persons for information

Sales			primediagnosics@wur.nl
Invoicing			primediagnosics.invoices@wur.nl
General information	José van Beckhoven	(+31) 317 480 603	jose.vanbeckhoven@wur.nl
Marketing and sales	Willem Stol	(+31) 317 480 561	willem.stol@wur.nl
	Iris Bergervoet	(+31) 317 488 610	iris.bergervoet@wur.nl
Packaging and shipping	Fianne Lubbers	(+31) 317 486 614	fianne.lubbers@wur.nl
Prime Support	José van Beckhoven	(+31) 317 480 603	jose.vanbeckhoven@wur.nl
	Sander Grapendaal	(+31) 317 484 462	sander.grapendaal@wur.nl

ELISA	Anke Clerxk	(+31) 317 480 663	anke.clerkx@wur.nl
	José van Beckhoven	(+31) 317 480 603	jose.vanbeckhoven@wur.nl
IIF	José van Beckhoven	(+31) 317 480 603	jose.vanbeckhoven@wur.nl
Inocula	Petra van Bekkum	(+31) 317 486 615	petra.vanbekkum@wur.nl
Lateral Flow Devices	Antje de Bruin	(+31) 317 480 631	antje.debruin@wur.nl
	Gabriel Lorencini Fiorin	(+31) 317 489 357	gabriel.lorencinifiorin@wur.nl
Luminex	Sander Grapendaal	(+31) 317 484 462	sander.grapendaal@wur.nl
	Henry van Raaij	(+31) 317 480 943	henry.vanraaij@wur.nl
TaqMan	Sander Grapendaal	(+31) 317 484 462	sander.grapendaal@wur.nl
	Henry van Raaij	(+31) 317 480 943	henry.vanraaij@wur.nl

Abbreviations

AP	alkaline phosphatase
CO	coating
DAS	double antibody sandwich
ELISA	enzyme linked immuno sorbent assay
IIF	indirect immunofluorescence
IIF	indirect immunofluorescence
LFD	later flow device
mL	milliliter

NC	negative control
OD ₄₀₅	optical density at 405 nm
PC	positive control
pv.	pathovar
subsp	subspecies
VPC	lyophilized positive control
mL	milliliter

Antisera Virus

The antisera listed below are available for ELISA and Luminex xMAP, LFD format on request. Ordering information and details of these formats can be found in the webshop of Prime Diagnostics: shop.wur.nl/primediagnosics.

Pathogen (virus)	Acronym
Alfalfa mosaic virus	AMV
Alstroemeria carlavirus	AICV
Alstroemeria flower-banding virus	AIFBV
Alstroemeria mosaic virus	AIMV
Alternanthera mosaic virus	AltMV
Andean potato latent virus/Andean potato mild mosaic virus	APLV/ APMMV
Andean potato mottle virus	APMV
Apple chlorotic leafspot virus	ACLSV
Apple stem grooving virus	ASGV
Arabis mosaic virus	ArMV
Bean common mosaic necrosis virus	BCMNV
Bean common mosaic virus	BCMV
Bean yellow mosaic virus	BYMV
Beet western yellows virus	BWYV
Bell pepper mottle virus	BePMV
Calibrachoa mottle virus	CbMV
Carnation etched ring virus	CERV
Carnation latent virus	CLV
Carnation mottle virus	CarMV
Carnation necrotic fleck virus	CNFV
Carnation ringspot virus	CRSV
Carnation vein mottle virus	CVMV
Cherry leafroll virus	CLRV
Chrysanthemum virus B	CVB
Cucumber green mottle mosaic virus	CGMMV
Cucumber mosaic virus	CMV
Cymbidium mosaic virus	CymMV
Dasheen mosaic virus	DsMV
Hosta virus X	HVX
Impatiens necrotic spot virus	INSV
Iris yellow spot virus	IYSV
Kalanchoe mosaic virus	KMV
Leek yellow stripe virus	LYSV
Lettuce big vein associated virus	LBVaV
Lettuce mosaic virus	LMV

Pathogen (virus)	Acronym
Melon necrotic spot virus	MNSV
Mirafiori lettuce big vein virus	MiLBVV
Odontoglossum ringspot virus	ORSV
Onion yellow dwarf virus	OYDV
Pea early-browning virus	PEBV
Pea seed-borne mosaic virus	PSbMV
Pelargonium flower-break virus	PFBV
Pelargonium line pattern virus	PLPV
Pepino mosaic virus	PepMV
Pepper mild mottle virus	PMMoV
Plum pox virus ('sharka')	PPV
Potato leafroll virus	PLRV
Potato virus A	PVA
Potato virus M	PVM
Potato virus S	PVS
Potato virus V	PVV
Potato virus X	PVX
Potato virus Y	PVY
Prune dwarf virus	PDV
Shallot yellow stripe virus	SYSV
Squash mosaic virus	SqMV
Strawberry latent ringspot virus	SLRSV
Streptocarpus flower-break virus	SFBV
Tobacco mild green mosaic virus	TMGMV
Tobacco mosaic virus	TMV
Tobacco ringspot virus	TRSV
Tomato aspermy virus	TAV
Tomato black ring virus	TBRV
Tomato brown rugose fruit virus	ToBRFV
Tomato bushy stunt virus	TBSV
Tomato mosaic virus	ToMV
Tomato ringspot virus	ToRSV
Tomato spotted wilt virus	TSWV
Tulip Virus X	TVX
Zucchini yellow mosaic virus	ZYMV

Antisera Bacteria and Oomycetes

The antisera listed below are available for IIF, other formats on request. Ordering information and details can be found in the webshop of Prime Diagnostics: shop.wur.nl/primediagnosics

Pathogen (bacteria)	Acronym
<i>Curtobacterium flaccumfaciens</i> pv. <i>flaccumfaciens</i>	Cff
<i>Curtobacterium flaccumfaciens</i> pv. <i>oortii</i>	Cfo
<i>Clavibacter insidiosus</i>	Ci
<i>Clavibacter michiganensis</i> subsp. <i>michiganensis</i>	Cmm
<i>Clavibacter michiganensis</i> subsp. <i>sepedonicus</i>	Cms
<i>Dickeya chrysanthemi</i> (<i>Erwinia chrysanthemi</i>)	Ech
<i>Erwinia amylovora</i>	Eam
<i>Pectobacterium atrosepticum</i>	Eca
<i>Pseudomonas chichorii</i>	Pc
<i>Pseudomonas syringae</i> pv. <i>lachrymans</i>	Psl
<i>Pseudomonas syringae</i> pv. <i>mors-prunorum</i>	Psm
<i>Pseudomonas savastanoi</i> pv. <i>phaseolicola</i>	Psph
<i>Pseudomonas syringae</i> pv. <i>pisi</i>	Pspi
<i>Pseudomonas syringae</i> pv. <i>porri</i>	Pspo
<i>Pseudomonas savastanoi</i> pv. <i>savastanoi</i>	Pssa
<i>Pseudomonas syringae</i> pv. <i>syringae</i> (strain specific)	Pssy
<i>Pseudomonas syringae</i> pv. <i>tomato</i>	Pst
<i>Rhodococcus fascians</i>	Rhf
<i>Ralstonia solanacearum</i>	Rsol
<i>Xanthomonas axonopodis</i> pv. <i>dieffenbachiae</i>	Xcd
<i>Xanthomonas phaseoli</i> pv. <i>phaseoli</i>	Xcph
<i>Xanthomonas citri</i> pv. <i>fuscans</i>	Xcphf
<i>Xanthomonas axonopodis</i> pv. <i>begoniae</i>	Xcb
<i>Xanthomonas campestris</i>	Xccam
<i>Xanthomonas axonopodis</i> pv. <i>vesicatoria</i>	Xcv
<i>Xanthomonas axonopodis</i> pv. <i>fragariae</i>	Xf
<i>Xanthomonas hyacinthi</i>	Xch
<i>Xanthomonas hortorum</i> pv. <i>pelargonii</i>	Xcp
<i>Xanthomonas arboricola</i> pv. <i>pruni</i>	Xpru


The antiserum listed below is available for LFD, other formats on request. Ordering information and details can be found in the webshop of Prime Diagnostics: shop.wur.nl/primediagnosics

Pathogen (oomycete)	Acronym
<i>Phytophthora</i>	Phyto

Special terms and conditions of Wageningen Plant Research (WPR) concerning Prime Diagnostics Products


- These "Special terms of delivery" of Wageningen Plant Research (WPR) and the "General Terms and Conditions" of Wageningen University & Research will apply to all offers and all agreements between WPR and client concerning Prime Diagnostics Products. In case of conflict between both terms, these "Special terms and conditions" takes precedence over the "General Terms and Conditions".
- Only webshop and e-mailed orders that have been accepted by WPR will be handled. Orders by phone are not handled.
- Complaints regarding deficiencies of the products delivered by WPR should be deposited within 2 months after the purchase and should be supported by relevant test results obtained with the standard protocol of Prime Diagnostics. In the event that WPR declares a complaint to be founded, WPR will be exclusively obliged to effect performance as agreed upon as yet or to refund the purchase price paid, at WPR's exclusive discretion.
- All agreements between WPR and the clients are governed by Dutch Law only. All disputes shall be handled exclusively by the competent court in The Hague, The Netherlands.
- All products are supplied under the condition that they are for the exclusive use by the client. They may not be sold, integrated into other commercial applications nor handed over to third parties without distribution or supply and license agreement.
- Antisera are biological products and differences in reactivity between batches may occur. Therefore any guaranty given by WPR is limited to the relative reactivity in a standard performed DAS-ELISA with a 10 times diluted positive control originating from Prime Diagnostics.
- By placing any order with WPR, the client declares that he/she has read and accepts these "Special terms of delivery" as well as the "General Terms and Conditions" of WPR and that he/she has complete knowledge of the current product information.
- Client shall use the material in appropriate containment conditions only for research purposes. In no event WPR shall be liable for any use by client of the material or any loss, claim, damage, or liability of whatever kind of nature, which may arise from or in connection with this agreement or the use, handling, storage or transport of the material. WPR shall be safeguarded by the recipient company against any claim regarding these matters. Any damage or loss to the material during transport is at the purchaser's full risk.
- In the event of any liability of WPR, this liability will be limited to the invoice amount for that part of the order to which the liability pertains, on the understanding that this amount shall in no event exceed the amount WPR in such case will receive from its liability insurance.
- Most of the conjugated antibodies are stabilized with bovine serum albumin (BSA). Due to European legislation (EU DIR 1774/2002 and 668/2004) the use of BSA is not allowed in fields of human or veterinary medicine, agriculture, food or cosmetics. By ordering antibodies from WPR client acknowledges that the ordered product, containing BSA, will be used exclusively for research and analytical purposes. If required a form for the written declaration can be obtained on request or downloaded from our website www.primediagnosics.com.
- Most products are usually available from stock and are shipped within two weeks on receipt of the order. However, to ensure timely delivery, orders should be placed 8 weeks prior to the desired delivery date. In the event of force majeure WPR will be entitled to suspend performance of the agreement or to terminate the agreement without recourse to the courts and without any liability towards the client. Force majeure on the part of WPR means any circumstance beyond the control of WPR, for example: strikes, fire, war, damage, transport difficulties, export obstructions, defaults on the part of suppliers and legal bars to manufacture or to supply the products.
- Payment should be made by wire transfer (Euro Base Payment applying code SHA (shared costs)) using IBAN of WPR and BIC/SWIFT of RABOBANK within 30 days after date of the invoice.
- Permits, import charges or client's formalities necessary for import by the client of the ordered products are not the responsibility of WPR. It is the responsibility of the client that the client takes care of this before placing an order. All necessary permits need to be in the possession of WPR before an order can be shipped and WPR cannot be held responsible for any damages or import charges that might occur from import problems. The client has to verify if the use of the products is allowed by its government.
- Prices are CPT (Incoterms 2020) in Euros (€), excluding VAT and all orders are charged with € 25.00 for handling. Orders exceeding € 5000.00 are free of handling costs.
- Specificity of the antisera is tested against strains known to be pathogenic at the time of the last actual testing. WPR cannot be held responsible for possible false positive or false negative results caused by newly emerging pathogens, developed resistance, pathogen strains or plant or matrix substances.
- The client has to verify the suitability of the products before purchasing the products. When the product information isn't clearly enough the client has to ask WPR for more information.
- For payments the following VAT number of Wageningen Plant Research have to be used:
 - Domestic : NL811383696B07
 - Foreign : NL806511618B01
- WPR does not accept the responsibility for any direct or indirect damage that might arise from the use of delivered products. WPR shall be safeguarded by the recipient company against any claim regarding the delivery and / or the use of the delivered products.

Schematic overview of the DAS-ELISA procedure



Schematic overview of the Luminex xMAP procedure

The bead mixtures used contain hundreds of beads of different bead addresses. Each bead address is conjugated to one (1) specific antibody for one (1) specific plant pathogenic organisms. This makes it possible to test for multiple pathogens simultaneously, each bead address - plant pathogen combination is specific for the pathogen it was designed for.



Sample	Plant pathogen				
	1 (bead 12)	2 (bead 23)	3 (bead 41)	4 (bead 71)	5 (bead 80)
1	2500	200	1356	230	28
2	27	30	29	33	31
n	x	x	x	x	x


Schematic overview of the TaqMan procedure

After DNA isolation specific primers and probe are added to the sample. The probe consists of quencher and a fluorescent molecule. When fluorophore and quencher are in proximity, quenching inhibits any fluorescence signals.


When the target is present the primers and probe will bind to the DNA (left) at a specified temperature. If the target isn't present the primers and/or probe will not bind (right).




Next, the temperature is decreased and polymerase will elongate the primers.




During elongation the polymerase will displace and cleave the probe, releasing a fluorescent molecule which will be measured (left). The fluorescence detected in the quantitative PCR is directly proportional to the fluorophore released and the amount of DNA template present in the PCR. If the probe did not bind, it is not cleaved and no fluorescent molecule will be released.



After increasing the temperature the DNA denatures and two single DNA strands are formed and the whole process of binding and elongation can start over again. Now with the doubled amount of targets.




The above steps can be repeated up to 40 times, and each repetition more fluorescent molecules will be released. In this way the signal is amplified exponentially. When plotted in a graph it will look as following:





Schematic overview of the Lateral Flow Device procedure

A sample of plant sap is applied to the strip. The control line confirms the test worked properly. A visible test line at the test line region indicates the presence of the pathogen. If the pathogen is not present in your sample, no test line will appear, indicating a negative result.

When pathogen is present in the sample



When pathogen is not present in the sample



Comparison between ELISA and Luminex xMAP

ELISA	Luminex xMAP
Workflow	
1. Coat microtitre plate	1. Sample extraction
2. Sample extraction	2. Add bead mixture to the wells of microtitre plate
3. Transfer samples to the wells of microtitre plate	3. Transfer samples to the wells of microtitre plate
4. Incubate and wash	4. Incubate and wash
5. Add secondary antibody-AP conjugate	5. Add secondary antibody mixture
6. Incubate and wash	6. Incubate and wash
7. Add reporter	7. Add streptavidin-R-phycoerythrin
8. Incubate and analyse	8. Incubate, wash and analyse
Multiplex	
No, multiplex not possible	Yes, up to 50 targets simultaneous detectable
Speed	
Depending on incubation times 16 to 24 hours	Finished in less than 4 hours
Consumables and labour input	
For each new target extra amounts of consumables and labour needed	Consumables and labour input not influenced by number of targets (pathogens)
Flexible	
No, coated or pre-coated microtitre plates needed	Yes, on demand the tests can be performed
Generally, whole plate has to be used	Depending on sample numbers partial plates can be used
Evaporation of samples in the wells at the edges of the plate	All wells of microtitre plate can be used
Reliability	
Duplos needed for reliable results	No duplos needed, for each targets more than 50 repetitions are measured

Accessing webshop and creating an account


To access the web site of Prime Diagnostics scan the QR code:



Or [click here](#) to go directly to our webshop

If you already created an account: login
If you are a first time user: create an account. You can follow the steps in the guide shown below.

Select "Create an Account" in the upper right corner of the webshop



This will open the following screen:

Create New Customer Account

Personal Information

First Name *

Last Name *

Allow remote shopping assistance ?

Sign-in Information

Email *

Password *

Password Strength: No Password


Confirm Password *

Create an Account

Fill in the required fields (marked with *), Next click "Create an Account".

After the account is created, please fill in the billing address under "My Account"

If the shipping address is different you need to provide this information as well.



- My Account
- My Orders
- My Downloadable Products

- Address Book
- Account Information
- Stored Payment Methods

My Account

Account Information

Contact Information

Iris Bergervoet
iris.bergervoet@wur.nl
[Edit](#) | [Change Password](#)

Address Book [Manage Addresses](#)

→ **Default Billing Address**

You have not set a default billing address.

→ **Edit Address**

Default Shipping Address

You have not set a default shipping address.

Edit Address

Add address details

- My Account
- My Orders
- My Downloadable Products

- Address Book**
- Account Information
- Stored Payment Methods

Add New Address

Contact Information

First Name *

Last Name *

Company *

Phone Number *

Address

Street Address *

VAT Number

Country *

State/Province

City *

Zip/Postal Code *

Click "Save address"
Repeat for shipping address if needed

[Save Address](#)

Once this is completed you will be informed by e-mail that you have created your account.

Before the first order can be made, and to confirm that the account was made by you, the account needs to be activated. Mail a reply to Prime Diagnostics (primediagnosics@wur.nl) to confirm that your account needs to be activated.

After activation you have full access to the webshop of Prime Diagnostics.

primediagnos^tics

Prime Diagnostics provides high quality reagents for the detection of plant pathogens in various formats:

- ELISA
- Immunofluorescence
- Inocula
- Lateral Flow Devices
- Luminex xMAP and xTAG
- TaqMan

Prime Support:

From starters to experienced users in the field of plant diagnostics we offer 'state of the art' training with emphasis on improvement of the overall quality of the existing laboratory.

- Companies who want to start up plant diagnostic laboratory activities, Prime Support, a service from Prime Diagnostics can assist in the first steps to start-up plant laboratory activities.
- To existing plant health diagnostics laboratory operations Prime Support offers training to improve skills in plant diagnostics.